


L i p o s o m e   A u t o M a k e r

# Automated Multi-functional Liposome Manufacturing Equipment

plus Programmable Bioreactor

Partnership Project Approved by the Ministry of Economy, Trade and Industry "New Coordination Support Project" on June 26, 2009

**One-touch operation, produces any type of liposomes automatically and repeatedly**



Reference: Liposome Engineering Laboratory, Inc.  
Hashimoto Electronic Industry Co., Ltd.

## Features and benefits of Automated Multi-functional Liposome Manufacturing Equipment

1. This is the only unique piece of equipment in the world that can prepare, with only vortex-mixing, not only MLV (multilamellar vesicles) but also LUV (large unilamellar vesicles) and GUV (giant unilamellar vesicles). SUV can be made automatically additional Sonication-Units.
2. With our development of an automated extruder that will be attached to this equipment, uniform LUV (for example, 100 nm size) will be able to be prepared and available for use in DDS (drug delivery system). This will result in easy preparation and development of liposomal pharmaceuticals for cancer therapy (cancer cells usually create angiogenic vessels with intercellular spaces around 100 nm between their epidermal cells).
3. This equipment can prepare large amounts of GUV, suggesting that this equipment is useful for developing artificial cells and the next generation types of liposome encapsulating viruses or micelles. In addition,, the prepared GUV can encapsulate very high amounts of low molecular drugs and siRNA as well as intact proteins such as antibody, enzyme, etc, which can also be used for DDS by reducing their size.


### Liposomes encapsulating pharmaceuticals, genes, nucleic acids, antibodies, enzymes, etc.

	Compound	Liposome	Concentration (mM)	Trapping efficiency (%)
Liposomes encapsulating pharmaceuticals, genes, nucleic acids, antibodies, enzymes, etc.	Calcein	GUV	50	60
		GUV→LUV	50	12
		LUV	50	12
	γ-Globulin	GUV	50	50
		GUV→LUV	50	24
		LUV	50	17
	Rh-dextran	GUV	50	32
		GUV→LUV	50	18
		LUV	50	17

This equipment reduces the time necessary for one cycle of preparation, makes multiple overnight-productions and has excellent reproducibility. Thus unskilled people can operate it and prepare liposomes easily although manual preparation is quite difficult.

## Functional Liposomes that can be Produced with the Automated Liposome Manufacturing Equipment


1. Liposomes
  - Multilamellar vesicles (MLV)
  - Small unilamellar vesicles (SUV)
  - Large unilamellar vesicles (LUV)
  - Giant unilamellar vesicles (GUV)
2. Liposomes encapsulating pharmaceuticals, genes, nucleic acids, antibodies, enzymes, etc.
3. Liposomes sensitive to temperature, pH, magnetism, ultrasound, etc.
4. Liposomes coated with PEG, saccharide chains, etc.
5. Liposomes with bound ligands, such as proteins, peptides, nucleic acids
6. Reconstituted liposomes, such as proteoliposomes and virosomes
7. Liposome vaccines




# Overview of Liposome Preparation Procedure

- Dissolve phospholipid in ether (lipid ether)
- Load the equipment with pertinent solutions
- Enter the amount and conditions of preparation for the liposome
- Press the preparation start icon on the touch panel

Automated preparation start



Working Model of Automated Multi-functional Liposome Manufacturing Equipment



- **Liposome vaccines encapsulating soluble antigens (MLV, SUV, and LUV)** : Inject aqueous solution-based soluble antigen.
- **Liposome vaccines encapsulating lipid-soluble antigens (MLV, SUV, and LUV)** : Inject organic solvent-based lipid-soluble antigen.
- **Reconstituted liposomes** : Pump to the test tube the liposome solution through the B2 port and the reactant solution through the B3 port, and then carry out vortexing under the appropriate temperature. (B2: Bottle2, B3: Bottle3 Reference: Page4 Equipment configuration)
- Amount of liposome to prepare: 2 to 20 mL per run by automated single-run preparation and by automated continuous-run preparation up to 10 consecutive runs




# Features of Automated Multi-functional Liposome Manufacturing Equipment


## Varied automated preparation programs

- Automated liposome preparation programs: Registration of 10 sets of preparation conditions each for up to 12 kinds is available.
- Bioreactor programs: Up to 15 user-created programs can be registered.
- Automated cleaning and drying programs: Details of operation, such as the cleaning liquid volume, drying duration, and number of runs, can be specified.

## Equipment functions and features



## Equipment configuration




# The Touch Panel Allows Simple Operation and Easy Creation of Programs

## Touch panel operations for MLV preparation


(1) Touch panel screen when power is turned on




(2) Liposome preparation is selected on the Process selection screen



(3) MLV is selected for the liposome to be prepared




(4) A desired program No. for MLV preparation is selected




Up to 10 sets of preparation conditions can be registered

(5) Enter the amount to be prepared, number of runs, vortex conditions, etc.




In a similar manner, enter the conditions of film preparation, MLV preparation vortex, cleaning, and drying, and then press the start button.

(6) The screen during MLV preparation



The MLV preparation program is started, displaying the screen during MLV preparation shown above. The running step is indicated by the lighting of a small circle at the beginning of the corresponding line.

(7) MLV is completed, and recovered in the bottle 1 (B1).



### Other items to be set

#### Film preparation:

Number of vortex revolutions, duration of the vacuum level 1, 2, 3, and 4, and temperature

#### MLV preparation:

Regular and reverse rotation durations, regular and reverse rotation interval durations, temperature, and total duration

#### Cleaning and drying:

Selection among cleaning and drying intensity levels L, M, and H



# Functions of Programmable Bioreactor

## Features

- Automated operation combining vortex mixer processing with evaporation functions and continuous ultrasonic homogenizer
- Automated execution of the process of mixing, dispersion, reaction acceleration, concentration, and dispensing of six different solutions under nitrogen atmosphere, which are defined in user-created programs
- Non-manual operation is supported for the program instructions by selecting icons on the screen, and 15 programs can be registered.
- Automated cleaning and drying using either purified water or ethanol
- A part of the program or the whole program can be repeated for a specified number of times.

## Features of component units

### •Ultrasonic homogenizer

A water-cooled continuous homogenizer adapter is incorporated, and six syringe pump units are used to deliver liquid for a wide range of applications, such as diffusion, emulsification, reaction acceleration, dispersion, and high-quality mixing and defoaming. Combination of continuous, batch, or circulation processing and vortex mixer processing is also available. Injection during vibration is also allowed.


### •Syringe pump

Stepping motors are used to control the syringe stroke length through 60 mm at the resolution of 24,000. Syringes are selectable from 1, 5, and 10 mL. Concurrent operation of syringe pumps at respective speeds is allowed, enabling synchronization of mixing timing and precision delivery of liquids. The syringe pumps can be used also as dispensers. Residual solutions remaining in the syringes are purged by nitrogen gas after syringe push-out.

### •Vortex mixer

Powerful mixing by 0 to 2,500 rpm of vortex revolutions and vibration of 9-mm amplitude are allowed, and test tube surface temperature is controlled by means of hot-air blowing and radiation thermometer. The vortex mixer can serve as a highly efficient evaporator by extending the surface area and evacuating the test tube. It can also serve as a continuous concentrator by combining syringe pumps.




## Operation instructions



### Basic operation instructions

Pertinent parameters, such as the source, the destination, and the amount of delivery, will be extracted for a selected basic instruction.

An example of created programs



## Specifications of Automated Multi-functional Liposome Manufacturing Equipment

### • Automated liposome preparation program (optional)

Kinds of liposome: SUV, MLV, LUV, GUV, Functional liposomes, Reconstituted liposomes, Liposome vaccines incorporating pathogens, Liposome vaccines encapsulating soluble antigens (MLV, SUV, and LUV), and Liposome vaccines incorporating lipid-soluble antigens (MLV, SUV, and LUV)  
Amount of preparation: Selectable among 2 to 20 mL per run  
Number of continuous runs: 1 to 10  
Preparation time: About 2 hours per run (including cleaning and drying process)

### • Bioreactor

Six syringe pump units are used to control the liquid delivery to the vortex mixer having the evaporator function and to the continuous ultrasonic homogenizer unit. It is a fully automated bioreactor that enables users to create as desired a series of preparation programs by selecting operation instructions for processing and liquid delivery available on the touch panel.

### • Operation

A touch panel is used for program selection and creation.

### • Configuration

Syringe pumps : Six units	Vortex mixer : 0 to 2,500 rpm (regular and reverse rotation) with evaporator and hot-air heating
Syringes : 1, 5, and 10 mL	Ultrasonic homogenizer unit : Output 200 W, and water-cooled continuous ultrasonic homogenizer adapter

### • Gas displacement

Automated operation completely free of air by means of nitrogen displacement of the whole system is supported.

### • Cleaning

Automated cleaning and drying using either purified water or ethanol (details of operation, such as the cleaning liquid volume, drying duration, and number of runs, can be specified)

### • Required utilities

Dry compressed air: 0.2 to 0.5 MPa, 20 L/min	Cooling water circulator
Nitrogen gas: 0.2 to 0.5 MPa or more. 25 L per run	Organic solvent recovery unit
Vacuum pump: 20 L/min, ultimate pressure 1kPa	

### • Dimensions


About 758 mm (W), 492 mm (D), and 690 mm (H) (excluding projections)

### • Mass

About 95 kg (excluding the ultrasonic homogenizer power supply)

### • Power source


100 VAC, 50/60 Hz, 15 A(115VAC, 220VAC)



## ■ Specifications of Utility Unit (Optional)

- Air compressor : Oil-less compressor, discharge pressure of 0.6 to 0.8 MPa, flow rate of 20 L/min
- Clean air system : Hollow fiber membrane air dryer plus filter regulator  
Atmospheric dew point of -10°C, nominal filtration rating of 0.01 µm,  
purge volume of 3 L/min
- Vacuum pump : Teflon diaphragm vacuum pump, pumping speed of 20 L/min,  
ultimate pressure of 1 kPa
- Organic-solvent recovery unit :  
Water-cooled solvent recovery
- Cooling water circulator :  
Temperature control range of -20 to 20°C (up to normal temperature)  
Circulation pumping speed of 12 L/min (max.)
- Control box : Continuous ultrasonic homogenizer adapter cooling-water control  
Hollow fiber membrane air dryer control  
Organic-solvent pathway control
- Dimensions : About 510 mm (W), 710 mm (D), and 1150 mm (H)
- Mass : About 50 kg
- Power source : 100 VAC, 50/60Hz, 10A (115VAC, 220VAC)

Note: A separate nitrogen gas cylinder is required for the nitrogen gas.



## Partnership Project Approved by the Ministry of Economy, Trade and Industry as a “New Coordination Support Project” on June 26, 2009

Project name: Manufacturing and marketing business of multi-functional equipment for automated manufacturing of liposomes, such as those used for research and development

<Core enterprise>

### Hashimoto Electronic Industry Co.,Ltd.

Research and development, manufacturing, and maintenance of the equipment

3866-12 Takasu-cho, Matsusaka City,  
Mie Prefecture, Japan 515-0104

Tel:0598-51-3111 Fax:0598-52-4300

E-mail: hdk@hashimoto-inc.co.jp

<http://www.hashimoto-inc.co.jp>



<Partner enterprise: Reference>

### Liposome Engineering Laboratory, Inc.

Research and development of application technologies of various types of liposome  
Manufacturing method studies, market development, and quality evaluation

1577 Kurimamachiya-cho, Tsu City,  
Mie Prefecture, Japan 514-8507

(in the university grounds of Mie University)

Tel:059-231-5326 Fax:059-231-5328

E-mail: [info@lel.co.jp](mailto:info@lel.co.jp) <http://www.lel.co.jp>



Research conducted in cooperation with Mie University